

Reporting practical session 29-11-2022



Student name: Vincent Belpaire
Student ID: 02009131

University of Ghent
Faculty of Engineering and Architecture
Bachelor of Science in the Biomedical Engineering

Abstract

A well made report of an experiment is crucial in research. It should contain all necessary information so if one would read it again a couple of years after the execution of the experiment one could redo it without many errors. In short, the report must allow repetition. This is only possible if all used parameters and methods are mentioned. And as a check, the results and conclusions should also be displayed.

This little article contains three reports. Each report relates to one of the practicals done at UZ Gent Campus: NTA, 3D cell culture and H&E.

1 NTA, Nanoparticle Tracking Analysis

1.1 aim

Calculating the size and concentration of particles in scatter and fluorescent mode for EV-characterisation.

1.2 material and methods

An LM10-series microscope, which utilizes the principles of light scattering and Brownian motion, was used to visualize the particles in the sample as highlighted points with highlighted tails. With the *NTA 3.4 - Sample Assistant Build 3.4.4 - SA* software the concentration and size distribution were calculated. Two types of samples were used. Both samples were rEVs but one had a scatter mode (SM) pre-treatment and one a fluorescent mode (FM) pre-treatment. Table 1 displays the environmental settings for both types of pre-treatments. All samples were made in a diluent of $\frac{1}{8000}$.

	scatter mode	fluorescent mode
Temperature/C	21.40 – 21.50	21.90 – 22.0
Viscosity/cP	0.966682 – 0.964370	0.95521 – 0.952941
Camera Type	sCMOS	sCMOS
Laser Type	Blue488	Blue488
Camera Level	13	16
Syringe Pump Speed/AU	20	20
Detection Threshold	3	3

Table 1: Environmental settings

For each treatment type three videos were recorded.

1.3 results and discussion

Only the concentrations in particles per ml are displayed. To have the absolute particle concentration, i.e. the absolute particle count, these concentrations should be multiplied with 8000 ml.

Table 2 shows the concentrations of different treatment-types. Using the averages of SM and FM to calculate the ratio of fluorescent particles gives $\frac{7.57E+08}{8.61E+08} = 0.879$. A high ratio indicates a well-executed separation in the EV-separation phase prior to the EV-characterisation which is the case in this experiment.

	video 1	video 2	video 3	mean	standard error
SM	8,14E+08	8,69E+08	9,00E+08	8,61E+08	2,53E+07
FM	7,51E+08	7,51E+08	7,45E+08	7,57E+08	9,17E+06

Table 2: concentrations (particles/ml)

Figures 1 and 2 displays respectively the size distribution of the different SM and FM samples.




Figure 1: Size distribution - scatter




Figure 2: Size distribution - fluorescent

1.4 conclusion

NTA allows the researcher to analyse a relatively high resolution of particles and with the use of fluorescence makes it possible to score the EV-seperation.

2 3D cell culture

2.1 aim

Analysing the size and morphology of 3D cell cultures for groups with high glucose (HG) and low glucose (LG) and comparing them statistically.

2.2 material and methods

The following steps were executed on a given data set of tiff-files.

ImageJ was used to determine the pixel to micrometer conversion which gave as result that $500\text{ }\mu\text{m}$ corresponds to 313 pixels, hence $1\text{ pixel} = \frac{500}{313}\mu\text{m}$.

With AnaSP data was extracted from the tiff-files. This data contains several parameters whereof two were chosen, *Circularity* and *FeretDiameterMax*, for a statistical compairinson between the two groups (LG and HG). These two parameters were chosen because the shape, as *Circularity*, is compared as well as the size, as *FeretDiameterMax*.

In the statistical analysis a significance level of 5% was used. For both parameters a independent t-test was performed to compaire if their means are equal, null hypothesis, or significantly different, alternative hypothesis.

Note that to use the independent t-test, homoscedasticity must be fulfilled.

2.3 results and discussion

In figure 5 (look at end of this article) the results of the independent t-tests are listed. For the *FeretDiameterMax* there is no significant dfference, hence the null hypothesis for this parameter was accepted. This

implies that the size of the cell can not be used to distinct between groups.

For the *Circularity* the condition of homoscedasticity is not fulfilled thus the results of the t-test can not be used. But if homoscedasticity is not fulfilled in the first place then must there be a significant difference between the two groups, hence *Circularity* is a good parameter to distinct both groups from each other.

2.4 conclusion

Some parameters are good for distinction between groups while others are not. With the appropriate software tools it is possible to determine which parameters are suitable for this role.

3 H&E, Hematoxylin/Eosin staining

3.1 aim

Analysing 2 slides from mice tissue for their structural composition by identifying their tissue type and occurring structures.

3.2 material and methods

The slides were prepared using the standard H&E method at the UZ Gent campus which will not be discussed in this report.

The name of the method refers to the two used stains, i.e. hematoxylin and eosin. Hematoxylin has a purple color and binds to the nucleus while eosin is rather pinkish and remain in the cytoplasm.

3.3 results and discussion

Figure 3 displays muscle tissue. Three regions are indicated as I, II and III and have the following observations:

I skeletal muscle

II fat tissue

III tumorous tissue




Figure 3: Muscle tissue

Figure 3 displays kidney tissue. The two circles, indicated with I, are glomeruli of the kidney.




Figure 4: Kidney tissue

3.4 conclusion

With H&E staining it is possible to visualize tissue structures so that it can be analysed by humans via a microscope. Besides the clear images presented, there is still need of expertise to have good identification of the different components. It is not always clear to identify the different tissue structures and if tissue of an ill-subject would be analysed it would get even more challenging.

4 SPSS output

Independent Samples Test										
		Levene's Test for Equality of Variances		t-test for Equality of Means						
		F	Sig.	t	df	Significance One-Sided p	Two-Sided p	Mean Difference	Std. Error Difference	95% Confidence Interval of the Difference Lower Upper
FeretDiameterMax	Equal variances assumed	1,646	,215	,423	19	,339	,677	10,082660869	23,852446244	-39,84108288 60,006404615
	Equal variances not assumed			,407	11,429	,346	,691	10,082660869	24,754990113	-44,15394216 64,319263902
Circularity	Equal variances assumed	16,396	<,001	6,106	19	<,001	<,001	,16502996386	,02702728407	,10846120818 ,22159871954
	Equal variances not assumed			5,815	9,252	<,001	<,001	,16502996386	,02837793375	,10109977282 ,22896015490

Figure 5: Independent t-test results